(EN) **Bacteria binding swabs**

Devic Description

Cutimed® Sorbact® Swab is a sterile, bacteria and fungi binding wound dressing. It consists of a green Sorbact® wound

contact layer, which allows passage of wound exudate into a secondary dressing. The dressing is folded into 8 layers. The dressing can be used in compression therapy.

Mode of Action

Cutimed® Sorbact® binds microorganisms, such as Staphylococcus aureus (including MRSA), Streptococcus species, Escherichia coli, Pseudomonas aeruginosa and Candida albicans, as shown in vitro. These microorganisms are removed from the wound each time the dressing is changed.

Intended Use

Cutimed® Sorbact® Swab is intended for use in the management of clean, contaminated, colonized or infected exuding wounds, such as surgical wounds, traumatic wounds, cavity wounds, fistulas, pressure ulcers, diabetic ulcers and foot and leg ulcers.

Contraindications

There are no known contraindications to the use of Cutimed® Sorbact® Swab. Attention is drawn to the warnings and precautions.

Warnings and Precautions

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The dressing should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with oily products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Cutimed® Sorbact® Swab is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross-contamination.

Do not re-sterilize.

Do not use if the peel pouch is already open or damaged, as the sterility of the dressing is guaranteed only when the package is unopened and undamaged prior to use.

Prior to radiation therapy, remove Cutimed® Sorbact® Swab. A new dressing can be applied following treatment.

Cutimed® Sorbact® Swab does not contain x-ray thread and should not be used in open surgery. Cutimed® Sorbact® Swab is magnetic resonance imaging (MRI) safe.

Instructions for Use

- Prepare the wound according to local clinical practice. If desired, gently press a Cutimed® Sorbact® Round Swab against the moist wound surface to remove superficial microbial contamination.
- Select an appropriate dressing size for the wound. The dressing can be unfolded to a suitable size and may overlap the wound margins if needed.
- Remove the dressing from the pouch using an aseptic technique.
- If cutting the dressing, use an aseptic technique. Discard any open and unused dressing.
- Direct wound contact is a prerequisite for the effective binding of microorganisms. Always apply the dressings in a way that ensures direct contact with the wound surface. Moist wound conditions are essential to ensure the efficacy of Cutimed® Sorbact® and to avoid potential sticking to the wound. In dry or low exuding wounds the use of Cutimed® Sorbact® gel dressings is recommended. Alternatively, Cutimed® Sorbact® Swabs can be combined with an amorphous hydrogel (e.g. Cutimed® Gel).
- In cavity wounds and fistulas, unfold and fluff up one or more dressings and line/fill the cavity or fistula. Avoid dense packing.
- Depending on the amount of wound exudate, Cutimed® Sorbact® Swabs can be combined with all types of hydroactive secondary dressings such as foams, superabsorbent dressings or alginates.
- If required, use a bandage, tape or a transparent film dressing for fixation.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days. The final decision on the therapy regimen needs to be taken by health care professionals.

Storage and Disposal

Keep dry and away from sunlight. Disposal should be made according to local environmental procedures.

Notice

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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Cutimed® é uma marca comercial registrada da BSN Medical GmbH.
Sorbact® é uma marca comercial registrada da ABIGO Medical AB.
Cutimed® je registrirana blagovna znamka družbe BSN medical GmbH.
Sorbact® je registrirana blagovna znamka družbe ABIGO Medical AB.
Cutimed® teregister trademark oleh BSN medical GmbH.
Sorbact® teregister trademark oleh ABIGO Medical AB.