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Cutimed® Epiona

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Instructions for use

Cutimed® Epiona

Wound dressing for local interactive wound treatment

Product Description

Cutimed® Epiona is a wound dressing for local interactive wound treatment. Cutimed® Epiona is a sterile, freeze-dried matrix consisting of collagen and calcium alginate. The native structured collagen fibers form a three-dimensional open-pore matrix. Collagen is obtained from bovine dermis and consists of collagen types I, II and V. Due to its good wet stability when moistened and its elasticity, Cutimed® Epiona can easily be applied to the wound. Cutimed® Epiona does not contain any chemical crosslinking agents.

Cutimed® Epiona is also distinguished by its hemostyptic properties. Cutimed® Epiona supports wound healing, is soft and easy to mould, and has a good absorptive capacity.

Indications

- Cutimed® Epiona is indicated for the treatment of all wounds healing by secondary intention which are free of necrotic tissue, e.g.:
 - Diabetic ulcers
 - Venous ulcers
 - Pressure ulcers
 - Ulcers caused by mixed vascular etiologies
 - Trauma wounds and surgical wounds

Mode of action

The open-pore structure of Cutimed® Epiona offers a high level of capillary action, making it particularly suitable for absorbing wound fluid. Along with the absorption of the wound exudate, necrotic material, bacteria and fibrin coatings being shed from the wound, inflammation-inducing proteases are captured, as are cytokines. This supports and accelerates formation of the granulation tissue. Migration of fibroblasts is induced and the synthesis of collagen is stimulated. Through the continuous supply of the defect with high-quality collagen during the granulation phase and at the start of the epithelisation phase, the neosynthesis of the body's own collagen fibrils and fibers is supported. By building up the wound bed, the necessary proliferation and migration of epidermal cells is supported. Cutimed® Epiona can thus be applied to a moderately bleeding wound surface.

Precaution

If there is a suspicion of a wound infection during treatment, a suitable antibacterial wound dressing should be used or a systemic treatment should be introduced.

Directions for use

Wound preparation

Prior to treatment with Cutimed® Epiona, dry and necrotic tissue has to be removed.

Note

After carrying out an enzymatic debridement, the wound should be cleaned carefully with a sterile physiological saline or Ringers' solution before Cutimed® Epiona is applied. To achieve an optimal effect, Cutimed® Epiona should be applied directly to the entire wound bed. In wounds with little or no exudate, Cutimed® Epiona should be moistened with a saline solution or Ringer's solution. The biologically degradable Cutimed® Epiona will be absorbed naturally by the body. Cutimed® Epiona Dressing must be covered with an appropriate secondary dressing in order to maintain a moist wound environment. After first application, the wound should be treated at intervals of up to 72 hours, depending on the amount of exudate. Cutimed® Epiona matrix which has not been absorbed by the body can remain in the wound.

Attention

Do not use if the individual package has been damaged/is open. For „Use by“ date, see printing on the packaging. Cutimed® Epiona must not be resterilized.

Contraindications

Cutimed® Epiona must not be used with a known intolerance or allergy to one or more components. There are no clinical data on the use of Cutimed® Epiona during pregnancy and lactation. Use of this product requires a careful risk/ benefit assessment by the treating physician and it should only be used if the advantages are clear.

Side effects

Despite the low antigenic potential of collagen which has been proven in the literature, an intolerance reaction cannot be fully excluded. Individual instances have been reported in the literature of pre-existing infections potentially being exacerbated.

Interactions with other agents

Cutimed® Epiona Dressing should not be used with chlorine-releasing antiseptics (e.g. chloramine), with protein-precipitating substances (tannic acid, silver nitrate) or caustic agents (iodine tinctures) that alter proteins. Disinfectants and tanning medicinal products should not be applied in the form of a moist dressing over the Cutimed® Epiona Dressing.

Shelf life and storage

Store Cutimed® Epiona dry at room temperature (not over 25° C). Cutimed® Epiona has a shelf life of 36 months. Do not use after expiry. The expiration date of this product is printed on the packaging. Cutimed® Epiona is sterile packed. In order to guarantee sterility, please make sure that the package is intact. Cutimed® Epiona from an already opened sterile package must be discarded. Do not resterilize Cutimed® Epiona as this can have an influence on the physico-chemical properties of the matrix

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Store out of reach of children!

Product and pack sizes		
[REF] 73229-00	5 cm x 5 cm	10 pieces
[REF] 73229-01	10 cm x 10 cm	10 pieces

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