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Product Description

Cutimed® Gelling Fiber is a soft, hydrophilic, sterile, nonwoven, gel forming pad or ribbon dressing composed of sodium carboxymethylcellulose and enhanced fibers for strengthening. Cutimed® Gelling Fiber is a conformable and highly absorbent dressing. Cutimed® Gelling Fiber absorbs wound fluid and transforms into a soft gel, which maintains a moist environment to support the bodies healing process. Cutimed® Gelling Fiber helps in the removal of nonviable tissue from the wound (autolytic debridement), without damaging newly formed tissue.

Intended purpose and indications

Under medical supervision Cutimed® Gelling Fiber is intended to cover dry to heavily exuding wounds and to absorb exudate.

- Cutimed® Gelling Fiber can be used for the local treatment of:
 - Leu ulcers, pressure ulcers (stage II-IV) and diabetic ulcers
- Surgical wounds (e.g. post-operative, wounds left to heal by secondary intent and donor sites)
- Partial thickness burns
- Exudate absorption in oncology wounds
- Traumatic wounds (e.g. abrasions and lacerations)

Material composition

Cutimed® Gelling Fiber is composed of sodium carboxymethyl cellulose and strengthening fibers made from high-density polyethylene and polyethylene terephthalate.

Directions for Use

Cutimed® Gelling Fiber is a moisture maintaining fiber dressing for dry to heavily exuding wounds. Product is cuttable to size and both sides of the product can be applied to the wound. Cutimed® Gelling Fiber can be left in place for up to 7 days, except in burns and donor sites where it can be left in place up to 14 days.

Please contact your healthcare professional if you are unsure whether the product is appropriate for you.

Preparation

Application and removal of the dressing should be in accordance with good clinical practice. Before applying the dressing disinfect hands. Prepare treatment environment and items for an aseptic dressing application (disposable gloves, sterile tweezers and scissors). If necessary, debride (e.g. with Cutimed® DebriClean) and irrigate the wound and remove leftover cleaning solution afterwards.

Application

Exuding wounds (excluding partial thickness burns)

- Choose the appropriate size of the dressing and cut to size, if necessary: Cutimed® Gelling Fiber should overlap the wound by 1 cm, as the dressing will shrink when it absorbs exudate. **(B1)**
- Apply Cutimed® Gelling Fiber to the wound and cover it with a moisture balancing dressing such as Cutimed® Siltec®, Cutimed® Sorbion® or other appropriate secondary dressing **(C)**.
- All wounds should be inspected regularly. Remove Cutimed® Gelling Fiber when clinically indicated (e.g. excessive bleeding, leakage or suspicion of infection).

Dry wounds

- Choose the appropriate size of the dressing and cut to size, if necessary: Cutimed® Gelling Fiber should overlap the wound by 1 cm, as the dressing will shrink when it absorbs fluid. **(B1)**
- Apply Cutimed® Gelling Fiber on the wound and wet with sterile water or saline over the wound area.
- Cover Cutimed® Gelling Fiber with a moisture retentive dressing such as Cutimed® Siltec® **(C)** to avoid drying out of the dressing and subsequent drying of the wound.
- All wounds should be inspected regularly. Remove Cutimed® Gelling Fiber when clinically indicated (i.e. leakage, excessive bleeding, or suspicion of infection).

Deep wounds and cavities

- Dry the surrounding skin thoroughly.
- Loosely pack Cutimed® Gelling Fiber ribbon or pad into the wound to allow swelling. When using Cutimed® Gelling Fiber ribbon, cut to appropriate length leaving a small overhang of 2 to 3 cm outside the wound for easy retrieval **(B2)**.
- Cover with a moisture retentive dressing such as Cutimed® Siltec® or other appropriate secondary dressing **(C)**.
- Partial thickness burns and donor sites**
- Choose the appropriate size of the dressing. For larger wound areas, Cutimed® Gelling Fiber dressings can be combined. Cutimed® Gelling Fiber should overlap the wound or other adjacent Cutimed® Gelling Fiber dressings by at least 5 cm, as the dressing will shrink when it is retracted **(B2)**.
- Cutimed® Gelling Fiber should be covered with sterile loosely woven gauze pad and appropriately secured **(C)**. In the post wound period (up to 4 days), high volumes of wound exudate may require frequent changes of Cutimed® Gelling Fiber.
- Change the cover regularly and inspect Cutimed® Gelling Fiber while it remains in place on the wound.
- Adherence to the wound bed of the Cutimed® Gelling Fiber is a desired characteristic. Non-adherence of areas of the dressing may indicate deepening of the wound or infection. Areas of the dressing may be cut away to enable inspection of the wound. The exposed areas should then be treated appropriately.
- As the wound re-epithelializes, Cutimed® Gelling Fiber will detach or be easily removed.

Removal (Dressing change)

- If removing Cutimed® Gelling Fiber is difficult, it should be fully saturated with sterile saline or water and removed slowly.
- Cleanse the wound after removal of Cutimed® Gelling Fiber to avoid residues with an appropriate sterile wound cleansing solution.
- Newly formed blood vessels may produce a blood-stained wound fluid.

Supportive measures

Appropriate supportive measures should be taken where indicated e.g. for venous leg ulcers use of graduated compression bandaging (e.g. JOBST® Ulcer Care)

Contraindications

- Cutimed® Gelling Fiber should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.
- Cutimed® Gelling Fiber is not intended for use as a surgical sponge.
- Cutimed® Gelling Fiber is not intended for use within internal body cavities or within closed wounds.
- Infection is not a contraindication to the use of Cutimed® Gelling Fiber.

Warnings and Precautions

Cutimed® Gelling Fiber should be used under medical supervision and if the treated condition deteriorates (e.g. development of infection), falls to improve or if a side effect is observed, consult a physician or an appropriate health care professional.

- Cutimed® Gelling Fiber should not be used if contaminated.
- Cutimed® Gelling is for single use only. Re-use may lead to increased risk of infection or cross contamination.

Size table

cm	inch
5 x 5 cm	2 x 2 in.
10 x 10 cm	4 x 4 in.
15 x 15 cm	6 x 6 in.
14 x 10 cm	3/2 x 4 in.
4 x 20 cm	3/2 x 8 in.
4 x 30 cm	3/2 x 12 in.
2 x 45 cm	3/4 x 18 in.
1 x 45 cm	2/5 x 18 in.

Storage and /or special handling conditions that apply

Store in cool dry place (10°C – 30°C/50°F – 86°F) and keep away from direct sunlight.

Disposal

Discard any unused portion of Cutimed® Gelling Fiber or used product according to local regulations, e.g. via regular household or clinical waste.

Notice

Any serious incident that has occurred in relation to this device should be reported to the BSN medical GmbH and the competent authority of your state.

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