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**Cutimed® Siltec® Sorbact®**

D

**Produktbeschreibung**

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ENG

**Product Description**

Cutimed® Siltec® Sorbact® is a sterile, bacteria-binding, absorbent polyurethane foam dressing for atraumatic dressing changes that contains super-absorbent stripes which absorb and lock additional wound exudate. The wound contact surface is bacteria-binding, constructed of Sorbact® mesh using the Sorbact® technology. The fixation border has a silicone coating. This minimises trauma and pain during dressing changes. The outer film is water-repellent, yet permeable to oxygen and vapour. The foam dressing and the super-absorbent stripes are designed to absorb and to lock away excess wound fluid inside the dressing, promoting a moist wound environment and minimising the risk of maceration. Cutimed® Siltec® Sorbact® is highly conformable and provides protection of the wound site. Cutimed® Siltec® Sorbact® leaves no residues in the wound. The dressing reliably retains wound exudate, even under compression.

**Indication**

Cutimed® Siltec® Sorbact® is indicated for the management of clean, colonized, contaminated or infected exuding wounds such as venous and arterial ulcers, pressure ulcers, diabetic gangrenes, surgical incisions, skin grafts and donor sites, lacerations or abrasions. Cutimed® Siltec® Sorbact® is recommended for wounds with moderate to high exudate levels.

**Contraindications**

No contraindications known.

**Warnings and precautions**

Cutimed® Siltec® Sorbact® is packaged for single use. Do not re-use or re-sterilize as there is a risk of transmission of body fluids or contaminated tissue between patients. Do not use if pouch is opened or damaged. Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide. Cutimed® Siltec® Sorbact® is not intended to be cut. Avoid overlapping two or more Cutimed® Siltec® Sorbact® dressings, always choose a dressing that is appropriate for the wound size. For external use only. Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms. The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing. Remove the dressing prior to radiation therapy. In wounds that show clinical signs of infection, a more frequent wound inspection is advised. In these cases, an appropriate systemic treatment should also be considered.

**Directions for use**

1. Clean the wound according to good clinical practice.
2. Make sure the surrounding skin is dry.
3. Choose an appropriate dressing size that covers the wound and overlaps the wound edges. (A)
4. Open the pouch and take out the dressing. Ensure sterility.
5. Remove protective film (B) and apply with the adherent side and the green wound contact layer towards the wound. (C)
6. Avoid stretching the dressing.
7. Depending on wound location it may be appropriate to apply additional dressing fixation.
8. Cutimed® Siltec® Sorbact® should be changed as often as necessary, dictated by the wound condition. As with all absorbent dressings, monitoring is required to ensure the dressing does not dry out and adhere to the wound. The dressing should be changed when saturation of wound fluid becomes visible.
9. Carefully remove any additional dressing fixation and gently remove Cutimed® Siltec® Sorbact® from the wound and discard.
10. If required, cleanse the wound prior to application of a new dressing.

**Storage**

Cutimed® Siltec® Sorbact® should be stored in dry conditions. Do not expose to direct sunlight.

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**BSN medical GmbH