## (ENG) Bacteria binding dressing pad

**Devce Description**

Cutimed® Sorbact® Dressing pad is a bacteria and fungi binding wound dressing. It consists of a green Sorbact® wound contact layer combined with an absorbent core and a white backing. Cutimed® Sorbact® Dressing pad absorbs and retains exudate, thereby reducing the risk of maceration and enabling a moist wound environment.

The dressing can be used in compression therapy.

**Mode of Action**

Cutimed® Sorbact® binds microorganisms, such as Staphylococcus aureus (including MRSA), Streptococcus species, Escherichia coli, Pseudomonas aeruginosa and Candida albicans, as shown in vitro. These microorganisms are removed from the wound each time the dressing is changed.

**Intended Use**

Cutimed® Sorbact® Dressing pad is intended for use in the management of clean, contaminated, colonized or infected wounds with moderate to high exudate levels, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic ulcers and foot and leg ulcers.

**Contraindications**

There are no known contraindications to the use of Cutimed® Sorbact® Dressing pad. Attention is drawn to the warnings and precautions.

**Warnings and Precautions**

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

Do not use on patients with a known hypersensitivity to any of the dressing components.

The green wound contact layer should always be applied in direct contact with the wound area to allow microorganisms to bind to the dressing.

Should the dressing adhere to the wound, moisten or soak the dressing to assist removal and avoid disruption of the healing wound.

Do not use in combination with oily products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Cutimed® Sorbact® Dressing pad is for single use only and should not be re-used. Re-use may lead to increased risk of infection or cross-contamination.

Do not re-sterilize.

Do not use if the peel pouch is already open or damaged as the sterility of the dressing is guaranteed only when the package is unopened and undamaged prior to use.

Do not cut the dressing.

Prior to radiation therapy, remove Cutimed® Sorbact® Dressing pad. A new dressing can be applied following treatment.

Cutimed® Sorbact® Dressing pad is magnetic resonance imaging (MRI) safe.

**Instructions for Use**

- Prepare the wound according to local clinical practice. If desired, gently press a Cutimed® Sorbact® Round Swab against the moist wound surface to remove superficial microbial contamination.
- Select an appropriate dressing size for the wound. The wound pad should overlap the wound margins by at least 2 cm.
- Remove the dressing from the pouch using an aseptic technique.
- Discard any open and unused dressing.
- Direct wound contact is a prerequisite for the effective binding of microorganisms. Ensure that the green contact layer comes into direct contact with the complete wound surface. Moist wound conditions are essential to ensure the efficacy of Cutimed® Sorbact® and to avoid potential sticking to the wound.
- If required, use a bandage or tape for fixation.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days. The final decision on the therapy regimen needs to be taken by health care professionals.

**Storage and Disposal**

Keep dry and away from sunlight. Dispose of waste according to local environmental procedures.

**Notice**

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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Cutimed® é uma marca registrada de BSN medical GmbH.
Sorbact® é uma marca registrada de ABIGO Medical AB.
Cutimed® je registrirana blagovna znamka družbe BSN medical GmbH.
Sorbact® je registrirana blagovna znamka družbe ABIGO Medical AB.
Cutimed® teregister trademark oleh BSN medical GmbH.
Sorbact® teregister trademark oleh ABIGO Medical AB.