Instructions for Use

Cutimed® Sorbion® Border

(EN)

Product description

Atraumatic self-adhesive, super-absorbent dressing. Cutimed® Sorbion® Border belongs to an assortment of strong exudate-absorbing wound dressings. Cutimed® Sorbion® Border is a sterile, hydroactive, gel-forming, self-adhesive wound dressing, offering atraumatic dressing changes. The dressing consists of a strong water-retaining superabsorbent core, a wound contact layer from a hypoallergenic polypropylene, a silicone adhesive border and a breathable, and a waterproof backing from polyurethane which allows to shower.

Characteristics

Cutimed® Sorbion® Border absorbs and reliably retains wound fluids within the superabsorbent core even under compression. This leads to a removal of wound healing inhibiting factors, including bacteria, proteases, and inflammatory mediators and minimizes the risk of leakage and maceration. At the same time, the dressing achieves a moist wound environment, which supports the body’s own wound healing. The high absorption capacity of Cutimed® Sorbion® Border enables prolonged wound rest compared to conventional wound dressings. Its enhanced fluid handling capacity reduces unpleasant odors. Nonviable tissue adheres to the dressing and are is gently removed at dressing change, supporting soft debridement activities. The top layer is especially flexible, waterproof and breathable, functioning as bacterial barrier. Its silicone adhesive border allows easy, secure and gentle fixation of the dressing. Cutimed® Sorbion® Border can be repositioned during application. The dressing can be atraumatically removed without leaving residues.

Intended Purpose

Cutimed® Sorbion® Border is intended for the treatment of exuding wounds with moderate to excessive levels healing by secondary intention such as decubitus, ulcus cruris, diabetic foot ulcers and other wounds healing through secondary intention, oozing laparotomy wounds, fistulas to the skin, exulcerating carcinomas, postoperative wound dehiscence and similar exuding wounds.

Please contact your healthcare professional if you are unsure whether the product is appropriate for you.

Contraindications

There are no known contraindications.

Warnings and Precautions

Do not use if the pouch is already open or damaged as the sterility of the dressing is guaranteed only when the pouch is unopened and undamaged prior to use. Cutimed® Sorbion® Border is for single use only. Re-use may lead to increased risk of infection or cross contamination. Do not re-sterilize. Cutimed® Sorbion® Border should be used under medical supervision. Please always follow the instructions of any product that is used in combination. Do not cut the dressing. It must remain intact and uncut. Before applying Cutimed® Sorbion® Border, make sure that the wound has sufficient amount of exudate. Avoid contact with the following areas: mucous membranes (risk of drying and adhesion) and eye area. If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate health care professional. A temporary increase in wound size may occur at the start of a normal healing process; consult a medical professional to determine the cause of any significant change in wound size, color, or odor. Cutimed® Sorbion® Border should not be used in combination with therapies involving increased heat development or in which high-energy radiation (UV, γ- or x-ray) is used, since this can compromise the mechanical integrity of the product.

Combination with other wound care products

Cutimed® Sorbion® Border should not be used in conjunction with petroleum or grease-based products, ointments or creams, as this may obstruct the pores of the dressing and reduce its efficacy.

Directions for Use

Application and removal of the dressing should be in accordance with good clinical practice.

Preparation:

- Before applying the dressing disinfect hands and use gloves. If necessary, irrigate the wound and remove leftover cleansing solution afterwards.
- Select a suitable dressing size, so that the dressing pad overlaps the wound edges by at least by 1cm/0.4inch, (Fig. A).
- Remove release liners, do not touch the wound pad (Fig. B).

- Apply the dressing towards the wound, (Fig. A). Apply the white, unprinted side of the product facing the wound.
- Gently press the adhesive border to the intact skin to secure the dressing, 6. Avoid stretching the dressing.

Additional dressing fixation:

Depending on wound location it may be appropriate to apply additional dressing fixation, e.g. to increase resistance to incontinence.

Combination with other products:

- Cutimed® Sorbion® Border may be used on infected wounds in addition to an antimicrobial therapy.
- If necessary, Cutimed® Sorbion® Border can be used as secondary dressing in combination with other dressings
 - with Cutimed® Sorbact® / Sorbact® in case of infected wounds,
 - with a wound contact layer on wounds with dry areas (e.g. Cuticel® Contact),
 - with a gel-based product on necrotic wounds (e.g. Cutimed® Gel).

Dressing change:

Cutimed® Sorbion® Border should be changed as often as clinically indicated, must be not fully saturated. Cutimed® Sorbion® Border can remain in place for heavily exuding wounds, up to 4 days. For wounds with moderate exudation, the dressing wearing time can be extended to 7 days.

Removal:

To remove Cutimed® Sorbion® Border lift one corner, gently remove the dressing and discard. Adherence of wound deposits to Cutimed® Sorbion® Border is a desirable effect, contributing to wound cleansing (soft debridement).

Storage and Disposal

Store dressing away from direct sunlight at ambient temperature and humidity. Discard used product according to local regulations, e.g. via regular household or clinical waste.

Notice

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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	Artwork ID: JBR16159.02	Review Loop: 4	Printing Colours:  Pantone 287 C  Cyan  Magenta  Yellow  Black
Project Title: C2503 Border IFU Rework		Date of Review Loop: 26.07.2021	
Contact: ArtworkManagementMedical@essity.com		Replaces Artwork ID: JBR16159.01	
Contact Address: BSN medical GmbH Schützenstraße 1–3 22761 Hamburg Germany		Technical Drawing for print ID: 465,5 x 286 mm	
Languages: EN, FR, DE, ES, IT, NL, PT, FI, SV, DA, NO, AR		Note Date of manufacture, LOT and expiry-date symbol will be part of variable coding.	

The following artwork was approved:

Asset Name	JBR16159.02
Component Name	JBR16159.02
Project Title	C2503 Border IFU Rework
Revision	2
Reason for Revision	Innovationsprojekt

Asset Status

Date (GMT)	Signed By	Role	Action	Comment
2021-06-30 09:25:26.0	Julia Schieback	Artwork Management	Artwork Submitted for Approval	I am the author of this document
2021-07-27 07:03:12.062	Helgo Schieber	Packaging Development	Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA)	I attest to the accuracy and integrity of this document
2021-07-26 14:39:41.948	Yvonne Link	Product Development	Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA)	I attest to the accuracy and integrity of this document
2021-07-26 13:37:01.738	Julia Otte	Regulatory Affairs Associate	Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA)	I attest to the accuracy and integrity of this document
2021-07-26 12:37:54.741	Ulrike Lehnigk	Global Medical Advisor	Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA)	I attest to the accuracy and integrity of this document
2021-07-27 08:10:34.0	Julia Schieback	Artwork Management	Final Sign Off	I attest to the accuracy and integrity of this document

THIS REPORT IS CREATED FROM ELECTRONIC SIGNATURE ACTIVITY WITH KALLIK AMS 360

CREATED ON 2021-07-27 08:10:34.0