# Cutimed® Sorbion® Sachet

Made and sterilized in Germany

[www.bsnmedical.com](http://www.bsnmedical.com)

Gebrauchsanleitung  
Instructions for use  
Mode d'emploi  
Instrucciones de uso  
Istruzioni per l'uso  
Gebruiksaanwijzing  
Manual de instruções  
Bruksanvisning  
Instruksjoner for bruk  
Brugsanvisning  
Käyttöohje

JBR11706.01

**Produktbeschreibung**

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Datum: Mai 2018

**Product Description**

The hydroactive Cutimed® Sorbion® Sachet wound dressing has an inner core with strongly water-retaining polymers (Superabsorber). The outer layer is made of ultrasonically sealed polypropylene fleece. The product does not contain any pharmacologically active substances and is free of binders or adhesives.

**Characteristics**

The inner core of Cutimed® Sorbion® Sachet consists of modified cellulose fibres and gel-forming super absorbant polymers (SAPs), which allows it to absorb large quantities of fluid. The SAPs and cellulose fibers are mechanically bound, so that upon contact with fluids Cutimed® Sorbion® Sachet maintain its stability. The dressing achieves a moist wound environment and at the same time good, retention of the absorbed fluid in the wound dressing. Improved control of bacterial colonisation of the wound is frequently seen and results in a reduction of unpleasant odours. Slough and toxins adhere to the dressing and are gently removed at dressing change, supporting soft debridement activities.

**Indications**

Formation of a sterile, mechanical barrier over the wound with absorption of exudates. Also applicable on exuding wounds such as decubitus, ulcus cruris, diabetic foot ulcers and other wounds healing through secondary intention, dehisced wounds, second degree burns, fistulae, drainages, ulcerating carcinomas and other types of exuding wounds.

Note: Cutimed® Sorbion® Sachet S Drainage, due to its specific shape, was developed for use with catheters, drainages, etc.

**Contraindications**

Avoid contact with the following areas: mucous membranes (risk of drying and adhesion), eye area, dry wound areas (particularly with simultaneous compression therapy), and areas with little to no exudate.

**Warnings/Precautions**

For single use only. Applying an used dressing on another patient can cause exchange of body fluids, leading to potentially harmful infections.

Do not resterilise. Do not reprocess. Do not use the product if the packing is damaged or has been inadvertently opened, unsterile dressings should be discarded.

Cutimed® Sorbion® Sachet is not designed, sold or intended for use except as indicated.

Only apply into tunnel-forming pockets such as decubitus ulcers under special monitoring. The product expands as it absorbs exudate and may adversely affect wound healing under certain circumstances.

Cutimed® Sorbion® Sachet products display high affinity to fluids. This may result in rapid fluid absorption from the wound site and adherence to the wound in drier areas. To avoid undesirable effects, users should apply a gel based product or a wound contact layer to the dry areas of the wound prior to applying Cutimed® Sorbion® Sachet. Medically desired fluids applied to the wound in the context of the treatment may also be absorbed. This should be taken into consideration when planning the treatment.

Cutimed® Sorbion® Sachet, depending on its size, can absorb several hundred grams of fluid. However, this can lead to a pressure load on the wound and the surrounding skin. Therefore, the dressing should be changed regularly at appropriate time intervals (see directions for use).

A temporary increase in wound size may occur at the start of a normal healing process; however, users should consult a medical professional to determine the cause of any significant change in wound size, color or odor.

**Application notes**

Cutimed® Sorbion® Sachet should only be used under the direction of a qualified healthcare professional.

Before applying the product, make sure that the wound is secreting a sufficient quantity of fluid. Cutimed® Sorbion® Sachet is recommended in wounds with moderate to high exudate levels. In less exuding areas, Cutimed® Sorbion® Sachet should be applied over an atraumatic wound contact layer.

During treatment, consider that viscous exudations are generally absorbed slowly and only to a limited extent.

1. Clean the wound following standard procedures. Please note the instructions below, regarding combination with other products.
2. Select a suitable dressing size to make sure the absorbent inner core overlaps the wound margins by approximately 1 cm / 0,4 inch. Do not open the polypropylene cover, nor cut the product to the right shape. Do not fold the device. Folding the product may affect maximum absorbability.
3. Depending on the condition of the wound surface, either apply Cutimed® Sorbion® Sachet directly as a primary dressing or over a wound contact layer. To guarantee the desired fluid absorption and reduce the risk of maceration in the wound environment, please ensure that there is sufficient contact between the wound surface and the inner core.
4. Secure the dressing after selecting an appropriate fixation for the condition of the wound surroundings and the skin. Films, medical tapes, bandages or any other suitable fixation dressing can be used. While fixing the dressing, take expansion during absorption into consideration. Cutimed® Sorbion® Sachet has been shown to be particularly effective in compression therapy.

5. Cutimed® Sorbion® Sachet should be examined regularly, particularly at the beginning of treatment. The dressing should be changed in a timely and regular manner in accordance with proper medical practice. Cutimed® Sorbion® Sachet should be replaced before becoming completely saturated and fluid penetrates the back of the dressing, at the latest after four days (seven days for Cutimed® Sorbion® Sachet Extra). This maximum wear time can be reduced, depending on the patient and the wound situation.

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

In some instances, wound deposits may adhere to the surface of the dressing, and reduce the product's absorbency. Such deposits are common and may be removed from the wound together with Cutimed® Sorbion® Sachet during dressing changes.

**Combination with other wound treatment products**

The manufacturer is not aware of any incompatibilities in combination with other wound treatment products. Nevertheless, combining the Cutimed® Sorbion® Sachet with products containing highly-oxidising ingredients (such as hydrogen peroxide or potassium permanganate), strongly alkaline ingredients or strongly acid ingredients should be avoided. Cutimed® Sorbion® Sachet should not be used in conjunction with petroleum or grease based products, ointments or creams, as this may obstruct the pores of the dressing and reduce its efficacy.

**Storage (and Disposal)**

Store in dry surroundings, keep away from sunlight and protect against high air humidity.

Keep out of the reach of children.

Disposal should be made according to local regulations.

Keep dry

Do not use if sterility is compromised or if the pouch is damaged

Do not resterilise

Do not reuse

Do not cut the product into pieces

Temperature limitation

Keep away from sunlight

Caution

**Sachet S**  
**Sachet XL**  
**Sachet Multi Star**  
**Sachet S Drainage**

Wounds with moderate to high levels of exudate

**Sachet Extra**

Wounds with high to excessive levels of exudate

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BSN medical GmbH