Cutimed® Sorbion® Sana

Made and sterilized in Germany

www.bsnmedical.com

Gebrauchsanleitung
Instructions for use
Mode d'emploi
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Manual de instruções
Bruksanvisning
Instruksjoner for bruk
Brugsanvisning
Käyttöohje

JBR11705.01

Produktbeschreibung

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BSN medical GmbH
Quickbornstrasse 24
20253 Hamburg - Germany
www.bsnmedical.com

Product Description

The hydroactive, atraumatic wound dressing Cutimed® Sorbion® Sana has an inner core with strongly water-retaining polymers (Superabsorber). This inner core is cocooned in two layers; a nonwoven interior layer of polypropylene fleece and an outer wound contact layer made of polyethylene film with three-dimensional perforations. Both layers are ultrasonically sealed. Cutimed® Sorbion® Sana does not contain any active pharmacological ingredients and is free of binders or adhesives.

Characteristics

The inner core of Cutimed® Sorbion® Sana consists of modified cellulose fibres and gel-forming super absorbant polymers (SAPs), which allows it to absorb large quantities of fluid. The SAPs and cellulose fibers are mechanically bound, so that upon contact with fluid Cutimed® Sorbion® Sana dressings maintain their stability. The dressing achieves a moist wound environment and at the same time good retention of the absorbed fluid in the wound dressing. Improved control of bacterial colonisation of the wound is frequently seen and results in a reduction of unpleasant odours. Slough and toxins adhere to the dressing and are gently removed at dressing change, supporting soft debridement activities.

Due to its 3-D structure and smooth surface, the cocooning polyethylene wound contact layer facilitates the fluid flow in the direction of the dressing, while offering few points of access for ingrowing tissue, consequently providing a protective layer against adherence to the wound bed. Cutimed® Sorbion® Sana is particularly suitable for wounds with a heterogeneous healing process.

Indications

Formation of a sterile, mechanical barrier to the wound with absorption of exudations. Cutimed® Sorbion® Sana is designed for use on chronic or traumatic exuding wounds, such as decubitus ulcers, ulcer cruris, diabetic foot ulcers and other and other wounds healing through secondary intention, dehisced wounds, as well as after burns, on donor sites and fragile skin, and epidermolysis bullosa (EB).

Contraindications

Avoid contact with the following areas: mucous membranes (risk of drying and adhesion) and eye area.

Warnings / Precautions

For single use only. Applying an used dressing on another patient can cause exchange of body fluids, leading to potentially harmful infections.

Do not resterilise. Do not reprocess.

Do not use the product if the packing is damaged or has been inadvertently opened, unsterile dressings should be discarded.

Cutimed® Sorbion® Sana is not designed, sold or intended for use except as indicated.

This product should not be used in combination with therapies involving increased heat development or in which high-energy radiation (UV, γ- or x-ray) is used, since this can compromise the mechanical integrity of the product.

Only apply to tunnel building pockets such as decubitus ulcers under special supervision. The product expands as it absorbs exudate and may adversely affect wound healing under certain circumstances.

Cutimed® Sorbion® Sana products display high affinity to fluids. This may result, if this product is not used as intended (e.g. prolonged use), in adherence to the wound or fusion with the surrounding tissue, despite its atraumatic properties. Medically desired fluids applied to the wound in the context of the treatment may also be absorbed. This should be taken into consideration when planning the treatment.

Cutimed® Sorbion® Sana, depending on its size, can absorb several hundred grams of fluid. However, this can lead to a pressure load on the wound and the surrounding skin. Therefore, the dressing should be changed regularly at appropriate time intervals (see directions for use).

A temporary increase in wound size may occur at the start of a normal healing process; however, users should consult a medical professional to determine the cause of any significant change in wound size, color or odor.

Application notes

Cutimed® Sorbion® Sana should only be used under the direction of a qualified healthcare professional.

Cutimed Sorbion Sana is recommended for the management of moderate to high exuding wounds including wounds with a heterogeneous healing process with low exuding areas. Despite the wound contact layer, Cutimed® Sorbion® Sana is not intended for the treatment of dry wounds.

During treatment, consider that viscous exudations are generally absorbed slowly and only to a limited extent.

1. Clean the wound following standard procedures. Please note the instructions below, regarding combination with other products.
2. Select a suitable dressing size to make sure the absorbent inner core overlaps the wound margins by approximately 1 cm / 0.4 inch. Do not open the inner nor outer polypropylene covers, nor cut the product to the right shape. Do not fold the device. Folding the product may affect maximum absorbability.
3. Depending on the condition of the wound surface, either apply Cutimed® Sorbion® Sana directly as the primary dressing or applied over other products. To guarantee the desired fluid absorption and reduce the risk of maceration in the wound environment,

please ensure that there is sufficient contact between the wound surface and the inner core.

4. Secure the dressing after selecting an appropriate fixation for the condition of the wound surroundings and the skin. Films, medical tapes, bandages or any other suitable fixation dressing can be used. While fixing the dressing, take expansion during absorption into consideration. Cutimed® Sorbion® Sana has been shown to be particularly effective in compression therapy.
5. Cutimed® Sorbion® Sana should be examined regularly, particularly at the beginning of treatment. The dressing should be changed in a timely and regular manner in accordance with proper medical practice. Cutimed® Sorbion® Sana should be replaced before becoming completely saturated and fluid penetrates the back of the dressing, at the latest of seven days. This maximum wear time can be reduced, depending on the patient and the wound situation.

If the treated condition deteriorates, fails to improve or a side effect is observed, consult a physician or an appropriate health care provider.

In some instances, wound deposits may adhere to the surface of the dressing, and reduce the product's absorbency. Such deposits are common and may be removed from the wound together with Cutimed® Sorbion® Sana during dressing changes.

Combination with other wound treatment products

The manufacturer is not aware of any incompatibilities in combination with other wound treatment products. Nevertheless, combining the Cutimed® Sorbion® Sana with products containing highly-oxidising ingredients (such as hydrogen peroxide or potassium permanganate), strongly alkaline ingredients or strongly acid ingredients/aliphatic or aromatic hydrocarbons, fats or halogenated hydrocarbons should be avoided. Cutimed® Sorbion® Sana should not be used in conjunction with petroleum or grease based products, ointments or creams, this may obstruct the pores of the dressing and reduce the its' efficacy.


Antiseptics based on Octenidin (up to 0.1%) or Polyhexanid (up to 2%) can be combined with the product. If other solutions are used for wound cleaning, they should be removed by rinsing with a sterile saline solution or one of the other solutions mentioned before applying Cutimed® Sorbion® Sana.

Storage (and Disposal)

Store in dry surroundings, keep away from sunlight and protect against high air humidity.

Keep out of the reach of children.

Disposal should be made according to local regulations.

	Keep dry		Wounds with moderate to high levels of exudate
	Do not use if sterility is compromised or if the pouch is damaged		superficial and deep wounds
	Do not resterilise		Do not reuse
	Do not cut the product into pieces		Temperature limitation 10°C / 50°F / 25°C / 77°F
	Keep away from sunlight		LATEX FREE FORMULA Not made with natural rubber latex
	Caution		STERILE EO Sterilised by Ethylene oxide
	essity		CE 0124
	® = registered trademark		BSN medical GmbH Quickbornstrasse 24 20253 Hamburg · Germany www.bsnmedical.com
			Date: May 2018

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The following artwork was approved:

Asset Name	JBR11705.01
Component Name	JBR11705.01
Project Title	C2008 - IFUs Cutimed Sorbion
Revision	1
Reason for Revision	new asset

Asset Status

Date (GMT)	Signed By	Role	Action	Comment
2018-05-30 09:30:33.0	Mirja Diehl	Artwork Management	Artwork Submitted for Approval	I am the author of this document
2018-06-20 10:19:01.657	Janina Burton	Regulatory Affairs HAM	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2018-06-04 13:32:40.474	Nils Blumenthal	Product Development	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2018-05-31 13:46:42.436	Helgo Schieber	Packaging Development	Artwork Approved : Electronic Signature	I attest to the accuracy and integrity of this document
2018-06-20 10:41:11.0	Mirja Diehl	Artwork Management	Final Sign Off	I attest to the accuracy and integrity of this document

THIS REPORT IS CREATED FROM ELECTRONIC SIGNATURE ACTIVITY WITH KALLIK AMS 360

CREATED ON 2018-06-20 10:41:11.0

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© = registered trademark
BSN medical GmbH Quickbornstrasse 24 20253 Hamburg - Germany www.bsnmedical.com